

RAD23B antibody



Catalog Number: 114514

Product name

RAD23B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAD23B Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer. Positive WB detected in Jurkat cells, HEK-293 cells, A431 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 58kd

Preparation

This antibody was obtained by immunization of RAD23B recombinant protein (Accession Number: NM_002874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114514(RAD23B antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114514(RAD23B antibody) at dilution of 1:50 (under 10x lens)