

SOX2 antibody

Catalog Number: 115514

Product name

SOX2 antibody

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

Antibody description

SOX2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 37-40kd

Preparation

This antibody was obtained by immunization of SOX2 recombinant protein (Accession Number: NM_003106). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115514(SOX2 antibody) at dilution of 1:600



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115514(SOX2 Antibody) at dilution of 1:5000