

Anti-RICTOR Mouse antibody



Catalog Number: 168514

Product name

Anti-RICTOR Mouse antibody

Specificity

Human, Monkey, Mouse

Antibody description

RICTOR Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human RICTOR expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot analysis using RICTOR mouse mAb against HeLa (1), PANC-1 (2), MOLT4 (3), and HepG2 (4) cell lysate.