

KATNA1 antibody



Catalog Number: 111914

Product name

KATNA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KATNA1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of KATNA1 recombinant protein (Accession Number: NM_001204076). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111914(KATNA1 antibody) at dilution of 1:500



IP Result of anti-KATNA1 (IP:Catalog No:111914, 4ug; Detection:Catalog No:111914 1:500) with mouse brain tissue lysate 3440ug.