

# RPS6KL1 antibody



Catalog Number: 114914

## Product name

RPS6KL1 antibody

## Specificity

Human, hamster; other species not tested.

## Antibody description

RPS6KL1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon tissue. Positive WB detected in CHO cells. Observed molecular weight by Western-blot: 59 kDa

## Preparation

This antibody was obtained by immunization of RPS6KL1 recombinant protein (Accession Number: BC004540). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



CHO cells were subjected to SDS PAGE followed by western blot with Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:50

# RPS6KL1 antibody



Catalog Number: 114914

---

(under 40x lens)