

# PPP2R3B Polyclonal Antibody



Catalog Number: 165914

## Product name

PPP2R3B Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to PPP2R3B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 356-575 of human PPP2R3B (NP\_037371.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions


WB 1:500 - 1:2000

## Validations



Western blot - PPP2R3B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PPP2R3B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Enhanced Kit .Exposure time: 30s.

Western blot - PPP2R3B Polyclonal Antibody

Western blot analysis of extracts of mouse skeletal muscle, using PPP2R3B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 30s.