

Anti-CLCN6 antibody



Catalog Number: 175914

Product name

Polyclonal

Anti-CLCN6 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat, Monkey

Applications

Antibody description

WB

Rabbit polyclonal antibody to CLCN6

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human clcn6. aa range:611-660

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:; Western Blot analysis of various cells using CLC-6 Polyclonal Antibody diluted at 1: 500