

LHCGR antibody

Catalog Number: 112214

Product name

LHCGR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LHCGR Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 100kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000233). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



WB result of Catalog No:112214(LHCGR



antibody).

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112214(LHCGR antibody) at dilution of 1:500