

# Profilin 1 antibody



Catalog Number: 114214

## Product name

Profilin 1 antibody

## Immunogen

[Human PFN1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Profilin 1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, Jurkat cells, NIH/3T3 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 15 kDa

## Preparation

This antibody was obtained by immunization of Profilin 1 recombinant protein (Accession Number: NM\_005022). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114214(PFN1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114214(PFN1 antibody) at dilution of 1:50 (under 10x lens)