

# Anti-PAEP antibody



Catalog Number: 175214

## Product name

Anti-PAEP antibody

## Specificity

Cow, Goat

## Antibody description

Rabbit polyclonal antibody to PAEP

## Preparation

This antigen of this antibody was beta-lactoglobulin

## Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

## Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB:1:500-2000

## Validations




Fig1: Sample: Goat milk at 30ug;; Primary: Anti-Beta-lactoglobulin at 1:300 dilution;; Secondary: HRP conjugated Goat Anti-Rabbit IgG(bs-0295G-HRP) at 1: 3000 dilution;; Predicted band size : 18kD; Observed band size : 18kD