

# Anti-Annexin A1 (5E4D8) Mouse antibody



Catalog Number: 168214

## Product name

Anti-Annexin A1 (5E4D8) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse, Rat, Monkey

## Applications

WB

## Antibody description

Annexin A1 (5E4D8) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human Annexin A1 beta protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Annexin A1 in A549, CHO-K1, COS7, C2C12, C6, 3T3 and HeLa cell lysates using Annexin A1 mouse mAb (1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.