

Anti-Annexin A1 (5E4D8) Mouse antibody



Catalog Number: 168214

Product name

Anti-Annexin A1 (5E4D8) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Annexin A1 (5E4D8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Annexin A1 beta protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Annexin A1 in A549, CHO-K1, COS7, C2C12, C6, 3T3 and HeLa cell lysates using Annexin A1 mouse mAb (1:1000 diluted). Predicted band size: 38 kDa. Observed band size: 38 kDa.