

# PBX1 antibody

Catalog Number: 113614

## Product name

PBX1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PBX1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in A2780 cells, HeLa cells. Positive IP detected in A2780 cells. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of PBX1 recombinant protein (Accession Number: NM\_002585). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:113614(PBX1 antibody) at dilution of 1:500



IP Result of anti-PBX1 (IP:Catalog No:113614, 4ug; Detection:Catalog No:113614 1:500) with A2780 cells lysate 1000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113614(PBX1 Antibody) at dilution of

# PBX1 antibody



Catalog Number: 113614

---

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)