

RBM42 antibody



Catalog Number: 114614

Product name

RBM42 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RBM42 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HEK-293 cells, HeLa cells, human liver tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 65-70kd

Preparation

This antibody was obtained by immunization of RBM42 recombinant protein (Accession Number: NM_024321). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114614(RBM42 antibody) at dilution of 1:500