

SSTR3 antibody

Catalog Number: 115614

Product name

SSTR3 antibody

Specificity

Human; other species not tested.

Antibody description

SSTR3 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Positive IHC detected in human pancreas tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001051). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115614(SSTR3 Antibody) at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:115614(SSTR3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:115614(SSTR3 Antibody) at dilution of 1:50 (under 40x lens)