

EEF1D antibody

Catalog Number: 110314

Product name

EEF1D antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EEF1D Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive IHC detected in human pancreas cancer tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of EEF1D recombinant protein (Accession Number: XM_047421433). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP


Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110314(EEF1D antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110314(EEF1D antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-EEF1D (IP:Catalog No:110314,

EEF1D antibody



Catalog Number: 110314

3ug; Detection:Catalog No:110314 1:300) with HeLa cells lysate 1200ug.