

# SLC22A18 antibody

Catalog Number: 115314

## Product name

SLC22A18 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SLC22A18 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, COLO 320 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 40-50kd

## Preparation

This antibody was obtained by immunization of SLC22A18 recombinant protein (Accession Number: XM\_047427033). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115314(SLC22A18 antibody) at dilution of



1:1000

IP Result of anti-SLC22A18 (IP:Catalog No:115314, 4ug; Detection:Catalog No:115314 1:600) with mouse brain tissue lysate 4000ug.