

TRIM38 antibody



Catalog Number: 116314

Product name

TRIM38 antibody

Specificity

Human; other species not tested.

Antibody description

TRIM38 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HEK-293 cells, human colon tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of TRIM38 recombinant protein (Accession Number: XM_024446303). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116314(TRIM38 antibody) at dilution of 1:800