

# CIDEC antibody



Catalog Number: 109314

---

## Product name

CIDEC antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CIDEC Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HeLa cells, L02 cells. Observed molecular weight by Western-blot: 29 kDa; 48 kDa

## Preparation

This antibody was obtained by immunization of CIDEC recombinant protein (Accession Number: XM\_047448685). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109314(CIDEC Antibody) at dilution of 1:1000