

Cynomolgus LTBR/TNFRSF3 (Fc Tag) recombinant protein



Catalog Number: 501714

General Information

Gene Name Synonym

Uncharacterized protein

Protein Construction

A DNA sequence encoding the cynomolgus LTBR (F6V995) (Met1-Met225) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 28

Molecular Mass

The recombinant cynomolgus LTBR is a disulfide-

linked homodimer. The reduced monomer comprises 439 amino acids and has a calculated molecular mass of 48.8 KDa. The apparent molecular mass of the protein is approximately 64-70 KDa respectively in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus LTBR / TNFRSF3 Protein (Fc Tag)
SDS-PAGE