

# ACSM2A Polyclonal Antibody



Catalog Number: 162714

## Product name

ACSM2A Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to ACSM2A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 418-577 of human ACSM2A (NP\_001295101.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:200 - 1:2000

## Validations



## Western blot - ACSM2A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACSM2A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.