

# Perforin antibody



Catalog Number: 113714

## Product name

Perforin antibody

## Immunogen

[Human Perforin Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Perforin Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse liver tissue. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of Perforin recombinant protein (Accession Number: NM\_001083116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113714(PRF1 antibody) at dilution of 1:500