

Anti-POLG Rabbit antibody



Catalog Number: 167714

Product name

Anti-POLG Rabbit antibody

Specificity

Human

Antibody description

POLG Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human POLG

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of POLG in HeLa,A549,HL-60,U2OS cell lysates using POLG Rabbit mAb(1:1000 diluted).Predicted band size:140kDa.Observed band size:140kDa.