

C2orf49 antibody

Catalog Number: 108714

Product name

C2orf49 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

C2orf49 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of C2orf49 recombinant protein (Accession Number: XM_047445805). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108714(C2orf49 antibody) at dilution of 1:500