

DRP1 antibody

Catalog Number: 110024

Product name

DRP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DRP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells, Transfected HEK-293 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 75-80 kDa

Preparation

This antibody was obtained by immunization of DRP1 recombinant protein (Accession Number: NM_012063). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110024(DNM1L,DLP1 Antibody) at dilution of



IP Result of anti-DNM1L,DLP1 (IP:Catalog No:110024, 4ug; Detection:Catalog No:110024 1:500) with mouse brain tissue lysate 3600ug.