

# TANK Polyclonal Antibody



Catalog Number: 165024

## Product name

TANK Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TANK

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human TANK (NP\_597841.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Western blot - TANK Polyclonal Antibody

Western blot analysis of extracts of HT-29 cells, using TANK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunofluorescence - TANK Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using TANK antibody. Blue: DAPI for nuclear staining.