

BHLHE41 Polyclonal Antibody



Catalog Number: 166024

Product name

BHLHE41 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to BHLHE41

Preparation

Antigen: A synthetic Peptide of human BHLHE41

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF, FC


Dilutions

WB 1:500 - 1:1000

IF 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - BHLHE41 Polyclonal Antibody

Western blot analysis of extracts of human brain, using BHLHE41 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.