

Anti-HORMAD1 antibody



Catalog Number: 176024

Product name

Anti-HORMAD1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to HORMAD1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human hormad1. aa range:321-370

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of various cells using HORMAD1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).