

Anti-BDNF Rabbit antibody

Catalog Number: 167024

Product name

Anti-BDNF Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BDNF Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human BDNF

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


Validations



Immunohistochemical analysis of paraffin-embedded human testis, using BDNF Antibody .



Western blot analysis of extracts of Human cerebellum lysate, using BDNF antibody.



Immunofluorescent analysis of HeLa cells, using BDNF Antibody.