

# BSA antibody



Catalog Number: 107024

---

## Product name

BSA antibody

## Specificity

Bovine; other species not tested.

## Antibody description

BSA Mouse Monoclonal antibody. Positive WB detected in Bovine serum tissue. Observed molecular weight by Western-blot: 67 kDa

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Bovine serum tissue were subjected to SDS PAGE followed by western blot with Catalog No:107024(BSA Antibody) at dilution of 1:4000