

Anti-TMEM85 Rabbit antibody



Catalog Number: 169024

Product name

Anti-TMEM85 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TMEM85 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human TMEM85

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/2000-1/10000

IP: 1/20

Validations



Western blot detection of TMEM85 in Jurkat,Rat Brain,C6,CHO-K1,HeLa cell lysates using TMEM85 Rabbit mAb(1:1000 diluted).Predicted band size:20kDa.Observed band size:20kDa.