

MAP4K2 antibody



Catalog Number: 112424

Product name

MAP4K2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MAP4K2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, Raji cells. Observed molecular weight by Western-blot: 85-91kd,58kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004579). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112424(MAP4K2 antibody) at dilution of 1:300