

# SMAGP antibody



Catalog Number: 115424

---

## Product name

SMAGP antibody

## Specificity

Human; other species not tested.

## Antibody description

SMAGP Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 25-28kd

## Preparation

This antibody was obtained by immunization of SMAGP recombinant protein (Accession Number: BC003379). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115424(Loc57228 antibody) at dilution of 1:600