

# Anti-INHBE antibody



Catalog Number: 176424

## Product name

Anti-INHBE antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to INHBE

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human inhbe. aa range:301-350

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1:: Western Blot analysis of L929 cells using Inhibin  $\beta$ -E Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100