

Anti-Bora Rabbit antibody

Catalog Number: 167424

Product name

Anti-Bora Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bora Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Bora

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of Bora in K562,C6 cell lysates using Bora Rabbit pAb(1:1000 diluted).Predicted band size:61kDa.Observed band size:61kDa.