

BCL2L12 antibody



Catalog Number: 108424

Product name

BCL2L12 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

BCL2L12 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive IHC detected in human breast cancer tissue. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 32 kDa




Preparation

This antibody was obtained by immunization of BCL2L12 recombinant protein (Accession Number: BC104005). Purification method: Antigen affinity purified.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:117097(BCL2L12 Antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:117097(BCL2L12 Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog

BCL2L12 antibody



Catalog Number: 108424

No:117097(BCL2L12 Antibody) at dilution of 1:50



Immunofluorescent analysis of (-20oc Ethanol) fixed PC-3 cells using Catalog No:117097(BCL2L12 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)