

Anti-MNX1 antibody



Catalog Number: 175824

Product name

Anti-MNX1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MNX1

Preparation

This antigen of this antibody was synthesized peptide derived from human protein . at aa range: 260-340

Formulation

Liquid, Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night.