

# DDX25 antibody

Catalog Number: 109824

## Product name

DDX25 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DDX25 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue, Jurkat cells. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 50-55 kDa

## Preparation

This antibody was obtained by immunization of DDX25 recombinant protein (Accession Number: NM\_001330438). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109824(DDX25 antibody) at dilution of 1:500



IP Result of anti-DDX25 (IP:Catalog No:109824, 4ug; Detection:Catalog No:109824 1:500) with mouse testis tissue lysate 4000ug.