

SEPX1 antibody

Catalog Number: 115124

Product name

SEPX1 antibody

Immunogen

[Human SEPX1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SEPX1 Rabbit Polyclonal antibody. Positive IHC detected in mouse kidney tissue.

Preparation

This antibody was obtained by immunization of SEPX1 recombinant protein (Accession Number: NM_016332). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded mouse kidney using Catalog No:115124(SEPX1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded mouse kidney using Catalog No:115124(SEPX1 antibody) at dilution of 1:50 (under 40x lens)