

TIGD1 antibody



Catalog Number: 116124

Product name

TIGD1 antibody

Specificity

human; other species not tested.

Antibody description

TIGD1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human brain tissue, human cerebellum tissue. Observed molecular weight by Western-blot: 68-70kd

Preparation

This antibody was obtained by immunization of TIGD1 recombinant protein (Accession Number: NM_145702). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116124(TIGD1 antibody) at dilution of 1:400