

# Anti-Calumenin (4C6) Mouse antibody



Catalog Number: 168124

## Product name

Anti-Calumenin (4C6) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human, Mouse

## Applications

WB

## Antibody description

Calumenin (4C6) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human Calumenin protein fragments expressed in E.coli.

## Validations

## Formulation


Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of Calumenin in SW480, mouse intestine and A431 cell lysates and using Calumenin mouse mAb (1:1000 diluted). Predicted band size: 37KDa. Observed band size: 45KDa.