

AP3M1 antibody

Catalog Number: 108124

Product name

AP3M1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP3M1 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Positive IP detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of AP3M1 recombinant protein (Accession Number: XM_024447939). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108124(AP3M1 antibody) at dilution of 1:400



IP Result of anti-AP3M1 (IP:Catalog No:108124, 4ug; Detection:Catalog No:108124 1:300) with NIH/3T3 cells lysate 1200ug.