

Anti-PDK2 Mouse antibody

Catalog Number: 169124

Product name

Anti-PDK2 Mouse antibody

Specificity

Human

Antibody description

PDK2 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PDK2 (AA: 178-404) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b


Applications

IHC-P, ICC/IF, FC


Validations



Flow cytometric analysis of HeLa cells using PDK2 mouse mAb (green) and negative control (red).





Western blot analysis using PDK2 mAb against human PDK2 (AA)



Immunohistochemical analysis of paraffin-embedded esophageal cancer tissues using PDK2

Anti-PDK2 Mouse antibody

Catalog Number: 169124



mouse mAb with DAB staining.

Western blot analysis using PDK2 mAb against HEK293 (1) and PDK2 (AA

Immunofluorescence analysis of HeLa cells using PDK2 mouse mAb (green). Blue