

# Anti-PDK2 Mouse antibody



Catalog Number: 169124

## Product name

Anti-PDK2 Mouse antibody

## Specificity

Human

## Antibody description

PDK2 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human PDK2 (AA: 178-404) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG2b


## Applications

IHC-P, ICC/IF, FC


## Validations



Flow cytometric analysis of HeLa cells using PDK2 mouse mAb (green) and negative control (red).



Western blot analysis using PDK2 mAb against human PDK2 (AA)




Immunohistochemical analysis of paraffin-embedded esophageal cancer tissues using PDK2

# Anti-PDK2 Mouse antibody




Catalog Number: 169124

---



mouse mAb with DAB staining.

Western blot analysis using PDK2 mAb against HEK293 (1) and PDK2 (AA



Immunofluorescence analysis of HeLa cells using PDK2 mouse mAb (green). Blue