

FARS2 antibody



Catalog Number: 110524

Product name

FARS2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FARS2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human cerebellum tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of FARS2 recombinant protein (Accession Number: XM_017010186). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110524(FARS2 antibody) at dilution of 1:500