

HOPX antibody

Catalog Number: 111524

Product name

HOPX antibody

Immunogen

[Human HOPX Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HOPX Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human placenta tissue, rat lung tissue. Positive IP detected in Mouse lung. Observed molecular weight by Western-blot: 8-12kd

Preparation

This antibody was obtained by immunization of HOPX recombinant protein (Accession Number: NM_139212). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



IP result of HOPX antibody (Catalog No:111524 for IP and Detection) with mouse lung tissue



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111524(HOPX Antibody) at dilution of 1:600