

# SP140L antibody



Catalog Number: 115524

## Product name

SP140L antibody

## Specificity

Human; other species not tested.

## Antibody description

SP140L Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 67-77 kDa

## Preparation

This antibody was obtained by immunization of SP140L recombinant protein (Accession Number: BC004921). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115524(SP140L Antibody) at dilution of 1:600



IP Result of anti-SP140L (IP:Catalog No:115524, 4ug; Detection:Catalog No:115524 1:500) with HeLa cells lysate 3200ug.