

Anti-AGTR2 antibody



Catalog Number: 176524

Product name

Anti-AGTR2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit monoclonal antibody to AGTR2

Preparation

This antigen of this antibody was synthetic peptide within c-terminal human angiotensin ii type 2 receptor.

Formulation

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IP

Dilutions

WB: 1:500-1:2,000

IP: Use at an assay dependent concentration.

Validations




Fig1.; Western blot analysis of Angiotensin II Type 2 Receptor on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control.; Lane 1: 293T cell lysate; Lane 2: Hela cell lysate