

Anti-ULK1 Rabbit antibody



Catalog Number: 167524

Product name

Anti-ULK1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ULK1 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human ULK1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


Validations



Western blot analysis of ULK1 expression in (1) HEK293 cell lysate; (2) PC12 cell lysate.



Immunohistochemical analysis of paraffin-embedded human heart, using ULK1 Antibody.



Immunofluorescent analysis of 293 cells, using ULK1 Antibody .