

CPSF1 antibody



Catalog Number: 109524

Product name

CPSF1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CPSF1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 160kd

Preparation

This antibody was obtained by immunization of CPSF1 recombinant protein (Accession Number: NM_013291). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109524(CPSF1 antibody) at dilution of 1:500